UNIVERSITY OF JAFFNA, SRI LANKA FACULTY OF ALLIED HEALTH SCIENCES

FIRST YEAR FIRST SEMESTER EXAMINATION IN BSc HONS (MLS)- 2020

MLSMT1132 MEDICAL LABORATORY TECHNOLOGY 1

Date: 14.02.2022 Time: 2 Hours

ANSWER ALL SIX QUESTIONS.

	0.00			
1.	1.1	Explain the working principle of an autoclave.	(20 Marks)	
	1.2	1.2 Briefly explain the causes leading to malfunctioning of an autoclave.		
			(30 Marks)	
	1.3	Describe the methods to confirm a successful sterilization process.	(20 Marks)	
	1.4	Explain the working principle of a laboratory digital pH meter and	give the steps	
		in calibration.	(30 Marks)	

- 2. 2.1 Explain the factors that affect the accuracy of an analytical balance. (30 Marks)
 - 2.2 Explain the applications and maintenance of the following equipment.
 - 2.2.1 Centrifuge (30 Marks)
 - 2.2.2 Water bath (20 Marks)
 - 2.2.3 Water distillation apparatus (20 Marks)
- 3. 3.1 Write short notes on
 - 3.1.1 Types of fire extinguishers and their respective applications. (30 Marks)
 - 3.1.2 Material Safety Data Sheet (MSDS). (20 Marks)
 - 3.2 3.2.1 Categorize the different biomedical wastes in a clinical laboratory.

(20 Marks)

3.2.2 Explain the methods which could be used to dispose the above-mentioned categories of wastes. (30 Marks)

- 4. 4.1 Name three blood borne viruses. (15 Marks)
 - 4.2 Describe the ways by which the above viruses can be transmitted in a clinical laboratory. (30 Marks)
 - Describe the preventive measures to be taken in a laboratory to prevent the blood borne virus infection. (35 Marks)
 - 4.4 List five (5) factors which can lead to infection in a person due to accidental exposure to blood borne viruses. (20 Marks)
- 5. 5.1 List the different classes of biosafety cabinet. (15 Marks)
 - 5.2 Give the working principle of each class of biosafety cabinet mentioned in Section 5.1. (30 Marks)
 - 5.3 Write two uses of each class of biosafety cabinet mentioned in Section 5.1.

 (25 Marks)
 - 5.4. What precautionary measures should an operator take before, while and after using the biosafety cabinet. (30 Marks)
- 6. 6.1 Describe how an image is formed in a phase-contrast microscope.

(30 Marks)

- 6.1.2 Give three applications of phase contrast microscope in a clinical laboratory. (20 Marks)
- 6.2 6.2.1 List the precautionary measures which should be taken to maintain a binocular light microscope. (20 Marks)
 - 6.2.2 Justify the precautionary measures mentioned in Section 6.2.1.

(30 Marks)

